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

# ADJUSTABLE DOUBLE FEMALE PYRAMID CONNECTOR 1K150

## Instructions for Use

### A PACK CONTENTS (Fig.1)

- a **1K150** double female pyramid connector including :
  - . a M36 threaded extension piece with 1.5 pitch **(1)** and its screw fittings
  - . a M36 tapped sleeve with 1.5 pitch **(2)** and its screw fittings
- 4 headless hexagonal socket screws M8 x 20 (not shown) [See Use §E => Fig. 5b]
- this leaflet of instructions for use.

### B USE AND CONNECTIONS (Fig.2)

- The **1K150** adjustable double connector is convenient for patients up to 100 kg, carried load included.
- The **1K150** connector should be installed between two male pyramid components **(3 and 4)**.

### C WHICH LENGTH SHOULD WE CUT OUT (Fig.3)

- Determine the build height of the product according to its location in the prosthesis (in the following operations this dimension will be called **Y**).
- Screw the threaded extension piece **(1)** into the tapped sleeve **(2)** until you reach the thread root. Measure dimension **X**.
- The length to be cut will be calculated as follows : **Length to be cut = X-Y**.

### D CUTTING TO LENGTH THE THREADED EXTENSION PIECE (Fig.4)

**Using the 5C413 cutting tool is highly recommended  
to cut the threaded extension piece to length.**

**Caution :** The cutting tool is 12 mm thick. This thickness indicates the minimal length of the re-cut threaded piece. NEVER re-cut the threaded piece so that it becomes shorter than 12 mm because this would lead the component to get fragile.

- Screw the threaded extension piece **(1)** into the **5C413** cutting tool so that the extension piece protrudes out side **F** of the tool by the length that should be cut.
- Install the assembly «threaded extension piece and cutting tool» into a vice and firmly tighten the jaws to lock the threaded extension piece.
- Use a hacksaw to cut the threaded extension piece **(1)** along side **F** of the **5C413** cutting tool.
- Loosen the jaws of the vice and remove the assembly «threaded extension piece and cutting tool» from the vice.
- Unscrew the threaded extension piece and remove it from the cutting tool.
- Carefully burr the threaded extension piece at the section level.

### E ASSEMBLING THE 1K150 CONNECTOR ONTO THE PROSTHESIS (Fig.5a / 5b)

- Re-screw the threaded extension piece **(1)** into the tapped sleeve **(2)** until you reach the thread root.
- Install the **1K150** connector onto the prosthesis.
- Check that the **(2a)** screws are in the sagittal plane and the **(2b)** screws in the frontal plane. If there are not so, you can slightly unscrew the tapped sleeve **(2)** so that the **(2a)** and **(2b)** screws be well positioned.






**Caution :** Never unscrew the tapped sleeve by more than a quarter turn. Otherwise the component would get fragile.

- Tighten the **(2c)** screw of the sleeve up to 7.4 Nm.
- Proceed to the alignment adjustments.
- If one of the screws **(1a)** or **(2a)** or **(2b)** is too short, replace it by one of the 20 mm long screws included in the pack (fig. 5b).
- Cover the pyramid locking screws **(1a / 2a / 2b)** with XC047 adhesive and tighten them up to a 15 Nm tightening torque.

### F RECYCLING

The **(1)** and **(2)** components are made of aluminum alloy and the **(1a)**, **(2a)**, **(2b)**, **(2c)** screws are made of steel. Each of these components should be recycled according to the laws in force.

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