Référence 7N51099

Français	Page 1
English	Page 32

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If an acronym belongs to the following list, it will be pronounced as indicated in the second column.

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ENGLISH VERSION
Reference 7N51099

Contents

1. Presentation	36
2. Technical characteristics	36
3. Use :	
- Keyboard	37
- Switch	39
- Memory	39
4. Setups : 4.1 Input :	40
KEYBOARD :	
- Type of keyboard.....	40
- Pressed time	41
- Released time	42
- Long time	42
- Repetition speed	42
- Beep of activation	42
- Sound volume of the beep	42
- Key action type	42
SWITCH :	
- Pressed time	42
- Released time	42
- Long time	42
- Repetition speed	42
- Beep of activation	42
- Sound volume of the beep	42
- Key action type	42
- Number of vertical scanning	42
- Number of horizontal scanning	42
- Type of scanning	42
- Example	43
PREDICTION :	
- Messages	44
4.2 Voice synthesis	45
4.3 Display	45
4.4 Miscellaneous	45
4.5 Version of software	46
4.6 Stop DIALO	46
5. Connections	47
6. Support.....	47
7. Problems ?.....	47 48

1. PRESENTATION

DIALO is a communication aid for the people who can't speak and who can use a keyboard or a switch.

DIALO is usable whatever the physical disability directly while pressing on the keys of the keyboard or while acting on an adapted switch.

The keyboard of DIALO is equipped with removable keyguard.

DIALO is equipped with two displays to check what is written and to an interlocutor to read the message (this second display can be stopped if necessary).

DIALO integrates a voice synthesis which will pronounce the message at the time wanted by the user.

The software of DIALO allows a dynamic prediction: each started word starts the proposal for 5. The prediction varies according to the using (the more one word is used, more quickly it is proposed)

DIALO can memorize until 10 000 characters of message than the user will place under codes as it wishes it. For exemple : "I am PIERRE and I communicate with this device" under the code : CD.

DIALO has a battery which allows an autonomous use from 4 to 12 hours according to the kind of use.

Finally DIALO is provided with a metal insert for fixing on a adaptable support.

! DIALO follows the europeans certifications ECM EN55011 and CE 1801.3. For that, it is provided with a transformer with a ferrite on the cable.

2. TECHNICAL CHARACTERISTICS

Sizes : 21 cms x 18cms x 8,5 cms

Weight : 1,250 KG

- Upper side :
- Dynamic display with a keyguard.
 - A display 4 rows of 40 characters.
 - A display 2 rows of 40 characters.
 - A speaker.
- Right side :
- A connexion for a transformer.
 - A connexion for a switch .
 - A connexion for maintenance(parallele port) .
- Left side :
- A serial port for maintenance.
- Lower side :
- A metal plate to install it on the adaptable support 7T21.

Transformer	9 Volts ; 1,33 Ampères
Size of keyboard	15 cm x 5,8 cm
Key size	12mm x 12 mm
Memory performance	10 000 characters for the prediction 10 000 characters for the messages

3. USE

! CAUTION :

During the first use, it is necessary to connect DIALO by the transformer delivered on the main power and to let it connected 24 hours in order to load the battery. The transformer has to be connected on the right side of DIALO.

! CAUTION : BE CAREFULL NOT TO USE ANOTHER TRANSFORMER UNDER PENALTY OF DETERIORATION OF DIALO AND OF RUPTURE OF GUARANTEE

When the transformer is connected, DIALO shows a small rectangle in bottom on the right of its display. The autonomy of the battery is planned for a use from 4 to 12 hours according to the kind of use . When the battery is low, DIALO in bottom on the right of the display "BAT » flashing. The load of the battery takes of 4 and 7 hours.

! CAUTION :


A written message cannot exceed 640 characters, in the event of going beyond, DIALO stops the writing and emits a double beep.

! THE VOICE SYNTHETIZER CAN SPEAK BY DIFFERENT WAYS WITH UPPER OR LOWER CASE. See Appendix.


3.1 Keyboard :

DIALO is supplied with a keyguard. If you wish to remove this keyguard it is possible to unscrew the two screws located under the product and which holds the plastic bar blocking the keyguard. If DIALO is switch off you have just to press any key to awake it. The standard keyboard of the country appears then: « azerty » for France, « qwerty » for many countries, « qwertz » for germany.

The standard keyboard is with lower case and the main functions keys.

Function key :

SPEECH : when you press the key , DIALO says what it is written on the display. this function is adjustable (see setups chapter).

MEMORY : this key is used to put and to find the messages in memory.

to

DIALO has a word prediction . when you begin a word , you have a proposal of 5 words , you can choose by using one of this keys .
This prediction is dynamic. Each new word is in memory automatically, more a word is used faster it is on the display

F6

You can have complete sentences in memory.

When you begin a word, DIALO will write on the forth line of the display the most used sentence beginning by this letter . you can confirm your choice by pressing this key F6.

•))

This key is to call somebody . The volume is adjustable (See setups).


With this key you can go to different keyboards : numeric , characters with accents.

CAUTION : This keys is for tempory action (when you choose for exemple a number you return after to the standard keyboard.)

Different kind of keyboards : press the key



First step :

	M	F1	F2	F3	F4	F5	F6			
•))	Q	Z	E	R	T	Y	U	I	O	P
	A	S	D	F	G	H	J	K	L	Z
HOME	←	└┘	X	C	V	B	N	M	↓	<u>123</u>



Second step :


	M	F1	F2	F3	F4	F5	F6			
.))	0	1	2	3	4	5	6	7	8	9
	,	.	é	è	ê	à	ù	ç	+	-
HOME	←	SPACE	←	→	/	SPACE	VOL-	VOL+	↓	<u>ABC</u>

Volume setup



New step :

	MEM	F1	F2	F3	F4	F5	F6			
.))	!	?	"	'	`	^	-	+	/	=
	&	()	,	.	:	;	"	%	€
HOME	←	└┘	@	*	<	←	→	>	↓	<u>123</u>

This key erases the displays.
This key can be (undo) to cancel the erase

Shift key for upper or lower characters.

123

Key to choose the numeric pad

3.2 Switch

All the functions described in the chapter KEYBOARD are accessible by switch by anybody who cannot physically use the keyboard. An action on the switch starts the scanning on the keyboard and a new action choose the right key. This way is adjustable (see Setups).

3.3 Memory

DIALO allows to put in memory any messages. For that, it is necessary to press twice on the key


then to press on the letters (or numbers or both) under which you wish to put in memory your message then space. It is preferable to choose codes of which it is easy to remember. For exemple : JP for "I'm. Joel PANSARD, I cannot speak but I understand you very well and I communicate with this device, thank you". WD for "I want to drink". WDH for "I want to drink a hot drink". WDC for "I want to drink a cold drink".

Once the setting is ok, it is possible to erase the display.

3 way to find these messages :

- 1) you have to press once the key, then on the code and space.
- 2) Press the code, then space (be careful to use words for the code)
- 3) Press the code with a long press on the last letter. See page 45


! CAUTION :

The length of the message cannot exceed 640 characters. In the event of going beyond DIALO blocks the writing and emits a double beep. In the event of full memory, DIALO writes in bottom, on the right of the display a M flashing.

ERASE A MESSAGE : press 2 times the key Memory then Backspace then the code to erase, then space.

4. SETUPS



DIALO gives access to very many adjustments described below ; some in direct access on the keyboard like the sound volume others in access by a code. The values of the adjustments factory are indicated between bracket.

Direct access to the following adjustments :

Sound Volume of the voice synthesis, by the keys VOL+ and VOL - on the Numeric keyboard : value of 0 ... 10 **(5)**

Nondirect access to the following adjustments :

The access to these adjustments is done by the combination of the key
then CALL, then 123.

All the setups are classified by menus for a faster selection. It is possible to validate or cancel a modification of setup and to leave the adjustments constantly.

4.1 Input

4.1a - Keyboard

- *Type of keyboard :*
 - alphabetical,
 - by frequency of the letters in language concerned for example EJARIN in French,
 - standard of the country or phonetics (***standard of the country***).

Alphabetic keyboard :

	MEM	F1	F2	F3	F4	F5	F6			
•))	a	b	c	d	e	f	g	h	i	j
	k	l	m	n	o	p	q	r	s	t
HOME	←	└┘	u	v	w	x	y	z	↓	<u>123</u>

By frequency :

	MEM	F1	F2	F3	F4	F5	F6			
•))	e	t	a	o	i	n	s	h	r	d
	l	u	c	m	f	g	y	p	w	b
└┘	←	HOME	v	k	x	j	q	z	↓	<u>123</u>

Phonetic :

	(M)	F1	F2	F3	F4	F5	F6			
•))	e	t	th	a	oy	o	ey	i	ow	n
	s	sh	h	r	d	j	l	u	k	ch
└┘	←	m	f	g	y	p	w	b	v	z

• *Pressed time :*

It's the adjustment of the delay of activation of the keyboard . This adjustment is very useful for the users who suffer from tremolos : value of 0,01s ... 1s per step of 0,1s then 1s ... 5s by step of 1s (0,1s).

- *Released time* : it's the time imposed between two input, very useful for the people reached of tremors. Value of 0,1s with 1s by step of 0,1s then 1s ... 5s by step of 1s **(0,1s)**.
- *Long time* : time required to have the automatic repetition of the keys.← ,→ and EFF: 0,5s ... 1s by step of 0,1s then 1s ... 7s by step of 1s **(0,5s)**.
- *Repetition speed* : 0,2s ... 1s by step of 0,1s **(0,5s)**.
- *Beep of activation* : YES or not **(YES)**.
- *Sound volume of the beep* : independent of the voice synthesis : 0 ... 10 **(5)**.
- *Key action type* : press or release (the long input are not then possible any more) this possibility is very significant for the users who cannot press on a key but which can make slip their finger on the keyboard after the guide finger was removed. in this case the action will be taken into account only when the finger leaves the keyboard.

4.1b - Switch :

- *Pressed time* :
it's the adjustment of the delay of activation of the keyboard. Tthis adjustment is very useful for the users who suffer from tremors : value of 0,1s ... 1s per step of 0,1s then 1s ... 5s by step of 1s **(0,1s)**.
 - *Released time* : it's the time imposed between two input, very useful for the people reached of tremors. Value of 0,1s with 1s by step of 0,1s then 1s ... 5s by step of 1s **(0,1s)**.
 - *Long time* : time required to have the automatic repetition of the keys.← ,→ and EFF: 0,5s ... 1s by step of 0,1s then 1s ... 7s by step of 1s **(0,5s)**.
 - *Repetition speed* : for the keys ←, →and EFF: 0,2s ... 1s by step of 0,1s **(0,5s)**.
 - *Beep of activation* : YES or not **(YES)**.
- Sound volume of the beep : independent of the voice synthesis: 0 ... 10 **(5)**.
- *Key action type* : Press or release (the long input are not then possible any more) this possibility is very significant for the users who cannot **press** on a key In this case the action will be taken into account only when the finger leaves the keyboard.
 - *Number of vertical scanning* : 0 ... 5.
 - *Number of horizontal scanning* : 1 ... 5 **(3)**.
 - *Type of scanning* : key by key or row/column (row/column).

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- *Scanning speed* : 0,1s ... 1s by step of 0,1s then 1s ... 5s by step of 1s (1s).

Automatic scanning: if the choice is YES the scanning starts again after each choice. YES or not (YES).

- *Auditory scan* : NO, BIP, KEYS (all pronounceable) (NO)

4.1c - Prediction :

- *Activation* : YES or not (YES).
- *Clear dictionnaire* : YES with request for confirmation or not (NOT)
- *Messages*:
 - Clear messages: YES with request for confirmation (NOT by defect) or not
 - See messages : you can see all the messages and their codes with the arrows

For example : code: HH Message: hello how are you?

4.1d – Message restitution :

- *Choice* : hold, space or standard. See page 40.

4.2 Voice synthesis

- Tone: 0 ... 10 **(5)**.
- Speed : 0 ... 10 **(5)**.
- Type: this choice makes it possible to decide when the voice synthesis speaks.
 With the factory setup DIALO speaks only when you press the key « speech »
 Character by character YES or not **(NOT)** ; Word by word (DIALO speaks after each word): YES or not **(NOT)** ; Sentence by sentence (DIALO speaks the punctuation '!', ';', '!' or '?'): YES or not **(NOT)**

4.3 Display

- **LCD user** : . Contrast: 0 ... 10 **(5)**
- **LCD interlocutor** :
 - Activation : YES or not **(YES)**, indeed it is possible to stop this display,
 some users want to use only the synthesis.
 - Contrast : 0 ... 10 **(5)**
- **Screen VGA** : it is the dynamic keyboard.
 Contrast : 0 ... 10 in visible beach **(5)**.
 TEXT Color : black,yellow,green, blu, red.
 DISPLAYColor : black,yellow,green, blu, red.
 light : de 1 à 10
 -
- **Backlight** : YES or not **(YES)**.

4.4 Miscellaneous

- **Call** :
 Type: choice among 6 sounds .
 Sound volume : independent of the voice synthesis: 0 ... 10 **(5)**.
- **Stand by time** :
 Time before DIALO stops after the last action on the keyboard or the switch
 : 1min ... 5min by step of 1min **(3)**.
- **Factory sets** : YES with request for confirmation (NOT by defect) or not.

4.5 Version of software

4.6 Stop DIALO

5. CONNECTIONS

Connection on the right side of DIALO :

Power connection 2,5mm internal and 5,5mm

Parallel port DB25 for maintenance

Connection jack 3,5mm mono female for switch, normally open

Connection on the left side of DIALO :

Connector RJ9 for serial port (maintenance)

6. SUPPORT

DIALO is equipped with a metal plate and can be fixed on our support 7T21.

7. PROBLEMS?

DIALO DOESN'T START WHEN YOU PRESS ON THE KEYBOARD :

Connect the transformer.

Check if it not a problem of contrast of the displays.

If the problem persists, check the power of the outlet then of the transformer: it must deliver 9 volts.

If the problem persists, please connect DIALO during one hour and again check.

If the problem persists, contact the maintenance service of your supplier.

DIALO STARTS BUT THE KEYBOARD DOESN'T WORK :

Connect the transformer.

If the problem persists press a longer time on the keyboard. if the keyboard reacts go into the setup menu and decrease the delay of activation of the keyboard.

If the problem persists, check if the keyguard is a the right place.

If necessary remove the keyboard to check.

If the problem persists contact the maintenance service of your supplier.

DIALO DOES NOT SPEAK :

Connect the transformer.

If the problem persists, check the volume of the voice with the keys VOL + or VOL -.

If the problem persists contact the SAV of your supplier.

THE SWITCH DOES NOT WORK :

To connect the transformer.

If the problem persists, press longer on the switch, if DIALO work please go to the setups to decrease the delay of activation of the switcher.

If the problem persists check the correct operation of the switch either on another device or with a tester.

If the problem persists contact your supplier.

THE PREDICTION DOES NOT WORK :

Check in the setups menus if it is switch on.

THE DISPLAY OF THE INTERLOCUTOR DOESN'T WORK

Check in the setups menus if it is switch on.